

Progesterone Receptor protein (STJ504712)

STJ504712

GENERAL INFORMATION

Product Type	Proteins
Short Description	Progesterone Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5241
Gene Symbol	PGR
Uniprot ID	PRGR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 320-370 on human Progesterone receptor protein.
Immunogen Region	320-370
Specificity	
Immunogen Sequence	