

## Histone H3 protein (Phospho) (S10) (STJ504717)

STJ504717

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Phospho-Histone H3-S10 protein was developed from and has a target region of S10. For use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">8350/8351/8352/8353/8354/8355/8356/8357/8358/8968</a>
<b>Gene Symbol</b>	<a href="#">H3C1.H3C2.H3C3.H3C4.H3C6.H3C7.H3C8.H3C10.H3C11.H3C12</a>
<b>Uniprot ID</b>	<a href="#">H31_HUMAN</a>
<b>Immunogen</b>	Phosphorylated synthetic peptide taken within amino acid region 1-50 on HH3 protein.
<b>Immunogen Region</b>	1-50
<b>Specificity</b>	
<b>Immunogen Sequence</b>	