

Phantom protein (STJ504721)

STJ504721

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phantom-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Mid region epitope synthetic peptide corresponding to unique amino acid sequences on D. melanogaster Phantom protein.
Immunogen Region	
Specificity	
Immunogen Sequence	