

PINK1 protein (STJ504728)
STJ504728**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PINK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65018
Gene Symbol	PINK1
Uniprot ID	PINK1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human PINK1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	