

PKC epsilon protein (STJ504730)

STJ504730

GENERAL INFORMATION

Product Type	Proteins
Short Description	PKC epsilon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5581
Gene Symbol	PRKCE
Uniprot ID	KPCE_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human Protein Kinase C epsilon type protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	