

cGKI protein (STJ504735)
STJ504735**GENERAL INFORMATION**

Product Type	Proteins
Short Description	cGKI-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5592
Gene Symbol	PRKG1
Uniprot ID	KGP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 609-659 on human protein kinase G alpha subunit.
Immunogen Region	609-659
Specificity	
Immunogen Sequence	