

Phospholipase A2 IIA protein (Phospho) (STJ504740)

STJ504740

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Phospholipase A2 IIA-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5320
Gene Symbol	PLA2G2A
Uniprot ID	PA2GA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-150 on human Phospholipase A2 membrane associated protein.
Immunogen Region	50-150
Specificity	
Immunogen Sequence	