

Phospholipase A2 protein (Phospho) (STJ504744)

STJ504744

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Phospholipase A2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9373
Gene Symbol	PLAA
Uniprot ID	PLAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-500 on human Phospholipase A2-activating protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	