

Lyn protein (Phospho) (STJ504745)

STJ504745

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Lyn-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4067
Gene Symbol	LYN
Uniprot ID	LYN_HUMAN
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 425-525 on human Lyn protein. Phospho Y491
Immunogen Region	425-525
Specificity	
Immunogen Sequence	