

MAP3K12 protein (Phospho) (STJ504749)

STJ504749

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-MAP3K12-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on MAP3K12 protein. Phosphorylated S298 & S302 residues.
Immunogen Region	270-320
Specificity	
Immunogen Sequence	