

MAPK14 protein (Phospho) (STJ504752)

STJ504752

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-MAPK14-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1432
Gene Symbol	MAPK14
Uniprot ID	MK14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	