

MERTK protein (Phospho) (STJ504753)

STJ504753

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-MERTK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65037
Gene Symbol	Mertk
Uniprot ID	MERTK_RAT
Immunogen	Synthetic peptide corresponding to amino acids 746-757 of Human MERTK. Sequence: CKKIYpSGDYpYpRQG
Immunogen Region	746-757
Specificity	
Immunogen Sequence	