

KIRREL protein (Phospho) (STJ504758)

STJ504758

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-KIRREL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170643
Gene Symbol	Kirrel1
Uniprot ID	KIRRR1_MOUSE
Immunogen	Phosphorylated synthetic peptide with phospo tyrosine 737 and 738 on mouse Kin of IRRE-like protein 1.
Immunogen Region	around pTyr737/738
Specificity	
Immunogen Sequence	