

PODXL1 protein (STJ504762)

STJ504762

GENERAL INFORMATION

Product Type	Proteins
Short Description	PODXL1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5420
Gene Symbol	PODXL
Uniprot ID	PODXL_HUMAN
Immunogen	Synthetic peptide corresponding amino acid region 200-250 on human PODXL protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	