

POMT1 protein (STJ504763)

STJ504763

GENERAL INFORMATION

Product Type	Proteins
Short Description	POMT1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10585
Gene Symbol	POMT1
Uniprot ID	POMT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-460 on human POMT1 protein.
Immunogen Region	400-460
Specificity	
Immunogen Sequence	