

PDE11A protein (Phospho) (STJ504773)

STJ504773

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE11A-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50940
Gene Symbol	PDE11A
Uniprot ID	PDE11_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	