

PDE11A protein (Phospho) (S117/S124) (STJ504774)

STJ504774

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE11A-S117/S124 protein was developed from and has a target region of S117/S124. For use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50940
Gene Symbol	PDE11A
Uniprot ID	PDE11_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on PDE11A protein. Phosphorylated Ser 117 & 124.
Immunogen Region	100-150
Specificity	
Immunogen	
Sequence	