

PDE3A protein (Phospho) (STJ504775)

STJ504775

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE3A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54611
Gene Symbol	Pde3a
Uniprot ID	PDE3A_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-300 on mouse PDE3A protein. Phosphorylated Serine 290 & 292 residues.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	