

PDE4A PKA site protein (Phospho) (STJ504780)

STJ504780

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE4A PKA site-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5144
Gene Symbol	PDE4D
Uniprot ID	PDE4D_HUMAN
Immunogen	Synthetic peptide corresponding to PKA phosphorylation site on PDE4A long form enzymes. The peptide was synthesized with phosphorylated serine residue.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	