

PDE4D PKA site protein (Phospho) (STJ504782)

STJ504782

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE4D PKA site-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5144
Gene Symbol	PDE4D
Uniprot ID	PDE4D_HUMAN
Immunogen	Phosphorylated synthetic peptide selective for PDE4D PKA site taken within amino acid region 50-100 on human PDE4D8 protein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	