

PDE4D ERK site protein (Phospho) (STJ504784)

STJ504784

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE4D ERK site-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5144
Gene Symbol	PDE4D
Uniprot ID	PDE4D_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on PDE4D4 variant gene.
Immunogen Region	
Specificity	
Immunogen Sequence	