

PGCA protein (Phospho) (STJ504788)

STJ504788

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PGCA-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4881
Gene Symbol	NPR1
Uniprot ID	ANPRA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Atrial natriuretic peptide receptor 1 protein.
Immunogen Region	500-550
Specificity	
Immunogen Sequence	