

PRP19 protein (Phospho) (STJ504794)

STJ504794

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PRP19-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246216
Gene Symbol	Prpf19
Uniprot ID	PRP19_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	