

PREX1 protein (STJ504797)

STJ504797

GENERAL INFORMATION

Product Type	Proteins
Short Description	PREX1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57580
Gene Symbol	PREX1
Uniprot ID	PREX1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 815-836 of Human PREX1. Peptide was covalently modified. Sequence: QEEDQADSAFPLLSLGPRLSLC
Immunogen Region	815-836
Specificity	
Immunogen Sequence	