

PREX2 protein (STJ504799)

STJ504799

GENERAL INFORMATION

Product Type	Proteins
Short Description	PREX2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80243
Gene Symbol	PREX2
Uniprot ID	PREX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 720-900 on PREX2 protein.
Immunogen Region	720-900
Specificity	
Immunogen Sequence	