

Prolactin protein (STJ504801)

STJ504801

GENERAL INFORMATION

Product Type	Proteins
Short Description	Prolactin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5617
Gene Symbol	PRL
Uniprot ID	PRL_HUMAN
Immunogen	N-epitope specific synthetic peptide taken within amino acid region 1-100 on human Prolactin protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	