

Prominin 1 protein (STJ504806)

STJ504806

GENERAL INFORMATION

Product Type	Proteins
Short Description	Prominin 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid regions 700-730 on CD133 protein
Immunogen Region	700-730
Specificity	
Immunogen Sequence	