

Src protein (Phospho) (STJ504822)

STJ504822

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Src-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6714
Gene Symbol	SRC
Uniprot ID	SRC_HUMAN
Immunogen	C-epitope phosphorylated synthetic peptide taken within amino acid region 400-500 on human Src protein. Phosphorylated Tyr-416 residue.
Immunogen Region	400-500
Specificity	
Immunogen Sequence	