

STAT3 protein (Phospho) (STJ504823)

STJ504823

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-STAT3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
Uniprot ID	STAT3_HUMAN
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 700-800 on human STAT3 protein.
Immunogen Region	700-800
Specificity	
Immunogen Sequence	