

PTIP protein (STJ504827)
STJ504827**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PTIP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22976
Gene Symbol	PAXIP1
Uniprot ID	PAXI1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1100 on human PTIP protein.
Immunogen Region	1000-1100
Specificity	
Immunogen Sequence	