

TORC2 protein (Phospho) (STJ504829)

STJ504829

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-TORC2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	200186
Gene Symbol	CRTC2
Uniprot ID	CRTC2_HUMAN
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequence on TORC2 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	