

TICAM2 protein (Phospho) (STJ504833)

STJ504833

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-TICAM2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
Uniprot ID	TCAM2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human TICAM2 protein. Phosphorylated serine at 16 position.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	