

TICAM2 protein (Phospho) (STJ504834)

STJ504834

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-TICAM2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
Uniprot ID	TCAM2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	