

Punt protein (STJ504839)

STJ504839

GENERAL INFORMATION

Product Type	Proteins
Short Description	Punt-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to C terminal amino acids 497-515 of Drosophila melanogaster PUNT protein. Peptides were amidated for conjugation. Sequence: QLLIKNHTNIDDAKESTNCSynthetic peptide corresponding to C terminal amino acids 497-515 of Droso
Immunogen Region	497-515
Specificity	
Immunogen Sequence	