

PODXL protein (STJ504842)

STJ504842

GENERAL INFORMATION

Product Type	Proteins
Short Description	PODXL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5420
Gene Symbol	PODXL
Uniprot ID	PODXL_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 510-560 on human PODXL protein.
Immunogen Region	510-560
Specificity	
Immunogen Sequence	