

Quiescin Q6 protein (STJ504843)

STJ504843

GENERAL INFORMATION

Product Type	Proteins
Short Description	Quiescin Q6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5768
Gene Symbol	QSOX1
Uniprot ID	QSOX1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human QSOX1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	