

Dexas1 protein (STJ504846)

STJ504846

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | Dexas1-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 19416 |
| Gene Symbol | Rasd1 |
| Uniprot ID | RASD1_MOUSE |
| Immunogen | Synthetic peptide: HSDLMYIREKTSVGSQAQDKERC corresponding to amino acids 256-278 of mouse Dexas 1 protein. |
| Immunogen Region | 256-278 |
| Specificity | |
| Immunogen Sequence | |