

RRF3 protein (STJ504850)
STJ504850**GENERAL INFORMATION**

Product Type	Proteins
Short Description	RRF3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1355-1405 on C. elegans RRF-3 protein.
Immunogen Region	1355-1405
Specificity	
Immunogen Sequence	