

RIC8A protein (STJ504854)
STJ504854**GENERAL INFORMATION**

Product Type	Proteins
Short Description	RIC8A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60626
Gene Symbol	RIC8A
Uniprot ID	RIC8A_HUMAN
Immunogen	Synthetic peptide taken within amino acids 500-535 on RIC8A protein
Immunogen Region	500-535
Specificity	
Immunogen Sequence	