

RICH2 protein (STJ504855)

STJ504855

GENERAL INFORMATION

Product Type	Proteins
Short Description	RICH2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	216831
Gene Symbol	Arhgap44
Uniprot ID	RHG44_MOUSE
Immunogen	Synthetic peptide corresponding to N terminal amino acids 1-60 of Mouse RICH2. Sequence: MKKQFNMRQ LANQTVGRA E KTEVLS E D LL QVEKRLE LVK QVSHST HKKL TACLQ GQQGA
Immunogen Region	1-60aa
Specificity	
Immunogen Sequence	