

RIPK1 protein (STJ504860)
STJ504860**GENERAL INFORMATION**

Product Type	Proteins
Short Description	RIPK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8737
Gene Symbol	RIPK1
Uniprot ID	RIPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 620-670 on human Receptor-interacting Serine/Threonine-protein kinase 1 protein.
Immunogen Region	620-670
Specificity	
Immunogen Sequence	