

RIPK2 protein (STJ504861)

STJ504861

GENERAL INFORMATION

Product Type	Proteins
Short Description	RIPK2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8767
Gene Symbol	RIPK2
Uniprot ID	RIPK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 410-460 on human RIPK2 protein.
Immunogen Region	410-460
Specificity	
Immunogen Sequence	