

RIPK3 protein (STJ504862)

STJ504862

GENERAL INFORMATION

Product Type	Proteins
Short Description	RIPK3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11035
Gene Symbol	RIPK3
Uniprot ID	RIPK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 461-511 on human Receptor-interacting Serine/Threonine-protein kinase 3 protein.
Immunogen Region	461-511
Specificity	
Immunogen Sequence	