

ROCK1 protein (STJ504867)

STJ504867

GENERAL INFORMATION

Product Type	Proteins
Short Description	ROCK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6093
Gene Symbol	ROCK1
Uniprot ID	ROCK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region spanning 850-900 on human ROCK1 protein.
Immunogen Region	850-900
Specificity	
Immunogen Sequence	