

## ROCK2 protein (STJ504868)

STJ504868

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	ROCK2-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	9475
<b>Gene Symbol</b>	ROCK2
<b>Uniprot ID</b>	ROCK2_HUMAN
<b>Immunogen</b>	Synthetic peptide spanning amino acid region 1300-1388 on human ROCK2 protein
<b>Immunogen Region</b>	1300-1388
<b>Specificity</b>	
<b>Immunogen Sequence</b>	