

ROR2 protein (STJ504870)

STJ504870

GENERAL INFORMATION

Product Type	Proteins
Short Description	ROR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4920
Gene Symbol	ROR2
Uniprot ID	ROR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-943 on human ROR2 protein.
Immunogen Region	900-943
Specificity	
Immunogen Sequence	