

S100 beta protein (STJ504876)

STJ504876

GENERAL INFORMATION

Product Type	Proteins
Short Description	S100 beta-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6286
Gene Symbol	S100P
Uniprot ID	S100P_HUMAN
Immunogen	Full length protein of human S100 beta/beta homodimer.
Immunogen Region	Full length
Specificity	
Immunogen Sequence	