

## EDG8 protein (STJ504879)

STJ504879

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	EDG8-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	60399
<b>Gene Symbol</b>	S1pr5
<b>Uniprot ID</b>	S1PR5_RAT
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 380-400 of rat S1P5R/EDG8 protein.
<b>Immunogen Region</b>	380-400
<b>Specificity</b>	
<b>Immunogen Sequence</b>	