

Sall4C protein (STJ504884)
STJ504884**GENERAL INFORMATION**

Product Type	Proteins
Short Description	Sall4C-protein if for use in research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human SALL4C protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	