

Saxophone protein (STJ504889)

STJ504889

GENERAL INFORMATION

Product Type	Proteins
Short Description	Saxophone-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids 107-125 of D. melanogaster SAX protein. Peptide was conjugated to KLH. Sequence: GCSTSPDQLPMICSQNSLK
Immunogen Region	107-125
Specificity	
Immunogen Sequence	